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The Philips equipment requires the following temperature and humidity conditions:

Requirement	Specification
Examination room	
Temperature	+18 to 22 °C
Maximum temperature change	5 °C/10 min
Relative humidity	40 - 70% (non-condensing)
Air conditioning capacity	2 kW
Control room	
Temperature	+18 to +24 °C (Human comfort) +10 to +35 (MR equipment)
Maximum temperature change	5 °C/10 min
Relative humidity	30 - 70% (non-condensing)
Air conditioning capacity	0.3 kW
Technical room	
Temperature	+15 to +24 °C
Maximum temperature change	5 °C/10 min
Relative humidity	30 - 70% (non-condensing)
Air conditioning capacity	8 kW (*)


(*) Heat load inside technical room is normally 2 kW. The heat load to air inside the technical room will increase dramatically if Air Cooled Cryo Cooler is activated. Therefore the hospital supplied air cooling must be able to deliver 8 kW cooling.


Note!
Heating, ventilation and air conditioning requirements concern all rooms (technical-, examination- and control room) and must be maintained 24 hours a day, 7 days a week.

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


The Philips equipment requires the following temperature and humidity conditions:


Requirement	Specification
Examination room	
Temperature	+18 to +24 °C
Maximum temperature change	4.1 °C/h
Relative humidity	35 - 70% (non-condensing)
Control room	
Temperature	+15 to +24 °C
Maximum temperature change	3 °C/hr
Relative humidity	35 - 70% (non-condensing)
Technical room	
Temperature	+15 to +28 °C
Maximum temperature change	5 °C/hr
Relative humidity	20 - 80% (non-condensing)

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

The Philips equipment requires the following temperature and humidity conditions:


Requirement	Specification
General requirements	
Temperature	+18 to +28 °C
Maximum temperature change	0,5 °C/minute.
Relative humidity	20 - 80%, non-condensing (non-condensing)
Examination room	
Heat emission	Average 1900 W
Control room	
Heat emission	Average 500 W
Technical room	
Heat emission	Average 2400 W

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
Heating, Ventilation & Air Conditioning (HVAC)

The Philips equipment requires the following temperature and humidity conditions:

Requirement	Specification
	Examination room
Temperature	+18 to +30 °C
Relative humidity	30 - 75% (non-condensing)
	Control room
Temperature	+5 to +35 °C
Relative humidity	20 - 75% (non-condensing)
	Technical room
Temperature	+10 to +40 °C
Relative humidity	30 - 75% (non-condensing)

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